



# PhysioWalk<sup>®</sup> 1.0

SAFE | COMPACT | EFFICIENT



# The Next Generation of Gait Trainers

## PhysioWalk® 1.0

PhysioWalk® 1.0 is a compact solution for effective gait training, designed to help maintain and improve walking ability and mobility. Thanks to its low step-up height, ergonomic access handrails, and integrated safety features, the system is ideally suited for use in home care, adult day care, and residential senior facilities.

The treadmill speed can be precisely adjusted even at very low speeds. The height-adjustable support bars, included as standard, can be optimally adapted to the user's body height.

## Key Features



**Low step-up height  
(8 cm, 11 cm with  
swing crossbars)**



**Touch panel  
for easy  
operation**



**Safety stop via dual  
light-barrier system**



**Large, high-visibility  
status display**



**Height-adjustable  
support bars**



**Speed safety mode  
for secure training**



*Gait trainer PhysioWalk® 1.0 in standard configuration*

## Optional Equipment


### **Forearm supports**

Swiveling and fully adjustable in height, width, rotation angle, and handgrip position to accommodate individual user needs.




### **Perturbation swing crossbars**

Designed to enhance activation of stabilizing musculature and improve sensorimotor function. The movement amplitude is infinitely adjustable.



# Technical Specifications PhysioWalk® 1.0



## Additional Models Available



PhysioWalk® 1.0 and 1.2 Basic



PhysioWalk® 1.2 Basic Side Panel



**PhysioWalk® 1.2**  
Walking surface 120 × 50 cm



**PhysioWalk® 1.8**  
Walking surface 180 × 50 cm

**FOR MORE INFORMATION**  
**TEL.: +49 75 62 – 9 72 60**  
**info@zebris.de, www.zebris.de**  
zebris Medical GmbH  
Am Galgenbühl 14, 88316 Isny, Germany

